

## Experiment to determine the concentration of a solution of hydrochloric acid

Each student will require the following apparatus and chemicals: —

Comboplate

Microtitration kit

Propette

Microspatula

White laminated sheet


Standard  $0.100 \text{ mol dm}^{-3}$  sodium carbonate solution

Approximately  $0.120 \text{ mol dm}^{-3}$  hydrochloric acid

Methyl orange indicator

The microburette apparatus used is shown below:

The base of the microstand should be positioned in well D2 and the two plastic arms placed in the way shown.



## Procedure

1) Rinse out your microburette with the solution.

hydrochloric acid

- a) Fill it up to the zero mark with the solution of hydrochloric acid.  
Make sure there are **no air bubbles** in the disposable tip.
- b) Place the microburette in the microburette stand as shown in the diagram opposite.
- c) Ensure that the disposable tip is not directly over the wells to be used.
- 2) Rinse out the 2 cm<sup>3</sup> pipette with the **standard** sodium carbonate solution.
- 3) Add **1 cm<sup>3</sup>** of sodium carbonate separately to wells F1, F2 and F3.
- 4) Add **1 drop** of methyl orange to each well.
- 5) Move the end of the microburette so that it is positioned over well F1.
- 6) a) Titrate the hydrochloric acid against the sodium carbonate solution by adding **a drop at a time**.  
b) After each addition **stir the mixture** with a microspatula.  
c) Continue adding the acid to the sodium carbonate until the **yellow** colour just changes to **orange**. If it goes red then you have gone past the end-point.  
d) Read the volume of solution from the microburette and record it in your results table on the other sheet.  
e) If necessary refill your microburette.
7. Repeat the titration with the sodium carbonate in wells F2 and F3 and keep on repeating until you get results which are within 0.010 cm<sup>3</sup> of each other\*.

**\* NOTE** as you get more skilled in your titration techniques try to get this figure down to 0.005 cm<sup>3</sup>

## Results

Volume of acid added	Titration 1	Titration 2	Titration 3
Final volume			

Starting volume			
Volume added			

Use the average of the two volumes added which are within  $0.01 \text{ cm}^3$  of each other.

Volume of acid added = \_\_\_\_\_  $\text{cm}^3$

**REMEMBER** the volume of the solutions must be quoted in  $\text{dm}^3$  and therefore your volumes have to be divided by 1000.  $1 \text{ cm}^3 = 1/1000 \text{ dm}^3 = 1 \times 10^{-3} \text{ dm}^3$

volume of acid added = \_\_\_\_\_  $\text{dm}^3$

The concentration of the sodium carbonate =  $0.05 \text{ mol.dm}^{-3}$

## Conclusions


The sodium carbonate solution is a **standard solution**. Explain what is meant by the term '**standard solution**'.

---



---

The equation for the reaction is: -



The equation shows that \_\_\_\_\_ mol. of HCl react with 1 mol. of sodium carbonate.

The number of moles of sodium carbonate (n sodium carbonate)

$$\begin{aligned}
 &= V_{\text{sodium carbonate}} \times C_{\text{sodium carbonate}} \\
 &= 0.1 \times 10^{-3} \times 0.05 \\
 &= 5.00 \times 10^{-5} \text{ mol.}
 \end{aligned}$$

$$n_{\text{HCl}} = \text{_____} \times 5.00 \times 10^{-5} \text{ mol.}$$

$$\text{Concentration of HCl} = n_{\text{HCl}} / V_{\text{HCl}}$$

=

The concentration of the hydrochloric acid

= \_\_\_\_\_ mol.dm<sup>-3</sup>